WillowWood[®]

ALPHA TrūSeal[™]

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Be sure to train the patient in the proper use and care of this product, and to provide the ***Alpha TruSeal Patient Instructions*** to the patient upon delivery of the product.

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Caution: Alpha Products can melt or burn if exposed to high temperatures or flame. Do not expose your liner product to these conditions. Doing so may result in ignition resulting in injury or death.

This Alpha product is intended for use on a single patient. Use of the product with multiple patients could lead to cross contamination between patients.

WARNING: To avoid danger of suffocation, keep this product and its packaging away from babies or children, and do not put the product on your head.

Do not expose the liner to items to which the patient is allergic. The liner could absorb the allergen and become a permanent source of exposure.

Silicone swells in contact with Thermoplastic Elastomer (TPE). Do not allow products made of TPE to contact the liner.

Consistent contact of the sealing rings against the socket wall is needed to maintain prosthetic suspension. If any dirt, lint, water, or other material disrupts the seal from contacting the socket wall, prosthetic suspension could become compromised.

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CONTENTS

SECTION	PAGE
Patient Criteria	4
Donning	5
Trimming	5
Shape Capture and Modification	6
Test Socket Evaluation	7
Definitive Socket Fabrication	7
Cleaning and Care.....	8
Troubleshooting	10
Warranty.....	12

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PATIENT CRITERIA

The liner will not be suitable for patient with any of the following conditions:

- Any patient whose limb circumferences do not fit within the size ranges for the liner.
- Any patient who is either not capable or not willing to comply with the cleaning and disinfecting procedures for the liner.
- Any patient who lacks the strength to roll the liner successfully up onto the residual limb.
- Any patient who lacks the ability to orient the liner properly on the limb.
- Any patient who has a skin condition which could be aggravated by covering the limb with a non-porous liner.
- Any patient with a residual limb shorter than 6 1/4 inches (160 mm), to ensure that all four sealing fins are inside the socket.
- Any patient experiencing substantial volume changes in their residual limb.
- Patients who are transtibial amputees.
- Patients with deep scarring or deep invaginations.

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DONNING

Do not apply any type of lotion or powder to the residual limb or to the liner, as these products could damage the liner. The residual limb should be clean, dry, and free of soap residue. Cover any open wounds or non-intact skin with a bandage or some other appropriate covering to prevent direct contact between the wound and the liner.

1. Invert the liner so that the fins are on the inside.
2. Place the end of the liner against the end of the limb. The orientation of the logo does not matter, but varying the orientation can prolong the life of the liner. (unless cut to fit trim lines)
3. Carefully roll the liner up onto the limb. To prevent irritation of the limb, make sure that there are no wrinkles or air pockets between the limb and the liner.

DO NOT PULL the liner onto the limb. Pulling the liner may stretch the skin and result in discomfort or skin irritation.

WARNING: If fabric socks or other wicking materials are used on the outside of the liner, do not cover the sealing fins. Covering the sealing fins will create a leak path and could cause loss of suspension.

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TRIMMING

Trim the liner to a shorter length if desired for patient comfort; however, keep in mind that maintaining a seal against the skin is critical for decreasing perspiration.

Use sharp scissors, a paper cutter, or a Hand-Held Cutting Wheel (WillowWood Part No. 700-TL002). If using the Hand-Held Cutting Wheel, place the liner on a cutting board or thick piece of plastic. Using a straight edge as a guide, roll the cutting wheel cross the liner.

Cut a smooth line. Avoid jagged scissor marks.

Using the medial socket trim line as the defining length, cut the liner either straight across or along the shape of the trimline.

Cutting the liner straight across allows the liner to be donned in any orientation. Cutting the liner to match the trimline is acceptable, as long as the patient is aware that the liner must be donned in the same orientation each time.

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SHAPE CAPTURE AND MODIFICATION

Fabricating a socket according to these guidelines is the most crucial element in prolonging the life of the Alpha Liner. If you do not follow the Total Surface Bearing concept, the life expectancy of the Alpha Liners will be significantly reduced and you will void the warranty.

Note: refer to the OMEGA® User's Guide and the OMEGA® Scanner 3D User's Guide for additional information.

Note: measure the residual limb length and circumference at 2" (50 mm) intervals.

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METHODS

Choose one of the following shape-capture methods according to patient criteria and residual limb assessment:

- OMEGA System
 - a. TF by Measurements
 - b. Scan and merge using existing or new templates.

Notes:

- Templates are included with the OMEGA software to assist in shape design and configuration. The software allows users to utilize the provided templates or add their own templates for reference.
- Because the residual limb tissue is held firmly in place by the liner, the patient may remain seated for the scanning process if desired.
- Traditional casting (using hand casting or a proximal brim system)

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MODIFICATION

Initially sand down the sealing rings until they are flush with the rest of the model. Reduce the shape in a gradient fashion from 4-5% at the proximal end to 0-1% distally at the seals.

The amount of reduction should be adjusted based on physical examinations and complementary measurements of the patient.

To facilitate diagnostic fitting and to reduce air entrapment, pay attention to these areas when modifying:

- Lateral wall shape
- Medial wall flare
- Adductor channel shape
- Distal compression

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TEST SOCKET EVALUATION

1. Fabricate a conventional diagnostic socket with an auto-expulsion valve.
2. Have the amputee don the liner and the test socket. The amputee may need to shift their weight a few times to achieve distal contact.

TIPS FOR ASSESSING THE STATIC DIAGNOSTIC FIT

- For lateral contact under load/weight bearing:
 - Assess amount of contact
 - Adduct socket if necessary to further evaluate contact
 - Distal contact should be achieved without overloading the distal femur, with minimal effort when donning, without elongating the tissues.
- Analyze total contact:
 - Look for signs of liner blanching
 - Observe proximal tissue containment
 - Assess overall comfort and support
- Confirm that all four of the sealing fins are contacting the socket.

CORRECT:
All four fins are contacting the socket.

INCORRECT:
All four fins are not contacting the socket.

- Mark the socket to indicate any areas of discomfort or improper fit, then modify the socket accordingly.
- Maintain Total Surface Weight Bearing.

DEFINITIVE SOCKET FABRICATION

Fabricate a definitive socket using a seamless flexible inner socket with an auto-expulsion valve and a rigid outer frame. Position the valve so that it is distal to the sealing fins on the liner.

The final socket should have smooth, flowing contours and no sharp edges. The inside of the socket must be perfectly smooth with no sharp edges.

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CLEANING AND CARE

Proper hygiene is extremely important when using the liner. Please review these care instructions very thoroughly with the amputee.

Daily Cleaning

Thoroughly clean the side of the liner without the fins. Use lukewarm tap water and a mild body soap that does not irritate the skin. Apply the soap with a clean, soft cloth or sponge.

Rinse **all** the soap residue from the liner with water.

Dry the liner with a clean, lint-free cloth.


Place the liner, **with the fins on the outside**, onto the Drying Stand* or onto a clean, dry, flat surface.

Clean the residual limb with soap and rinse well.

**Available separately in most countries*

Notes:

Do not scrub the liner. Scrubbing can roughen the surface of the liner, which can then irritate the skin.

Caution: Drying the liner with the fins facing in will damage the liner.

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Weekly Disinfection

Place a small amount of ethyl or isopropyl alcohol on a soft, clean cloth.

Lightly wipe the side of the liner without the fins with the cloth for two minutes.

Rinse off any excess alcohol with water.

Dry the liner with a clean, lint-free cloth.

Place the liner, **with the fins on the outside**, onto the Drying Stand* or onto a clean, dry, flat surface.

**Available separately in most countries*

Notes: Do not immerse the liner in ethyl or isopropyl alcohol. Extended contact with large amounts of ethyl or isopropyl alcohol will stiffen the liner.

Daily Inspection of the Liner

Explain to the amputee that every liner should be inspected after each day of use for the following conditions:

- Any signs of damage or unusual wear, which may indicate changes in the prosthetic fit.
- Torn fins.
- Debris in the fins.

If any of these conditions are observed, the amputee should stop wearing the liner and contact you immediately.

Perspiration

Any amputee that engages in an activity that could result in excessive perspiration should be instructed to periodically remove the liner and dry off the liner and residual limb to prevent slippage of the liner on the limb, and to dry the inside of the socket.

Storage

When the liner is not being worn, it should be stored with the fins on the outside in a cool, dry place on the Drying Stand or in the liner's original box. If the liner will not be used for an extended period, cover the dry liner with a clean plastic bag.

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TROUBLESHOOTING

Problem	Probable Cause	Solution
Residual limb and liner rotate inside the socket.	Residual limb is extremely cylindrical, and socket is loose.	Tighten socket fit.
Top edge of liner is rolling down the leg.	Liner is too short.	Change to a longer liner.
	Liner is too tight at the proximal end.	Confirm that the amputee is wearing the correct size of liner.
Open wounds and non-intact skin.	Could be caused by a number of issues, including socket fit and specific patient variables.	Address the specific cause of the situation, and apply a bandage or other appropriate covering to prevent direct contact between the wound and the liner.
Amputee who has been wearing an Alpha Liner for a long time without sweating has suddenly started sweating again.	Air is getting between the liner and the residual limb, possibly due to residual limb shrinkage.	Check the fit of the liner and socket. Tighten up the socket or change to a liner of a different size.
	Liner has a hole in it.	Replace the liner.
Amputee who is a new Alpha Liner user is perspiring excessively	Body has not yet adjusted to the closed environment of the liner.	<p>The perspiration should subside after several weeks of Alpha Liner use.</p> <p>Any patient that engages in an activity that could result in excessive perspiration should be instructed to periodically remove the liner and dry off the liner and residual limb to prevent slippage of the liner on the limb.</p> <p>Some amputees have found the moderate application of unscented alum-based anti-perspirant to the residual limb prior to donning the liner to be an effective method of controlling excessive perspiration. Always follow label directions, and never use anti-perspirant if there are open sores on the residual limb. Note: non-alum-based products may damage the liner.</p>

(continued on next page)

Problem	Probable Cause	Solution
Skin irritation all over the residual limb.	Poor hygiene.	Review proper cleaning, rinsing, and disinfecting procedures with the amputee.
	Use of an irritating lotion, cream, powder, or soap.	Check to see whether the amputee has recently changed to a different cleaning product, or has started using any lotions, creams, or powders containing hydrocarbon oils or animal fats or oils. (For some amputees, antibacterial soap has caused skin irritations.)
	The liner or the socket is loose, possibly due to residual limb shrinkage	Tighten up the socket or change to a liner of a different size, style, or thickness to eliminate the looseness. To confirm that this is not an allergic reaction, perform a “patch test”: apply a piece of Alpha material somewhere else on the body and look for a reaction.
Skin irritation along the proximal edge of the liner.	Amputee is pulling the liner up onto the skin instead of rolling it.	Refer to proper donning instructions on page 5. Cutting the proximal edge in a wave pattern may also reduce the irritation.

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WARRANTY

Important Note: For the safety of you and your patients, WillowWood will not resell any returned Alpha Liners that may possibly have been worn.

Material Safety Data Sheet

A Material Safety Data Sheet (MSDS) for the Alpha TruSeal Liner is available from WillowWood upon request.

Warranty

WillowWood warrants each Alpha Liner from the date of fitting for a period of 12 months against defects in material and workmanship. There is no warranty, expressed or implied, for damage caused by abuse, mishandling or accident. The warranty is void if the fitting and routine care instructions are not followed. For your protection and ours, please thoroughly clean all used products prior to return.

Warranty Disclaimer

WillowWood warrants that each product manufactured will, at the time of delivery, be of workmanlike quality and substantially free of defects. WillowWood makes no other warranty, implied or expressed, and makes no warranty of merchantability or fitness for a particular purpose. This warranty shall terminate immediately upon an action to combine our products with other materials or in any manner change the nature of our products. The remedy is repair or replacement or credit for the failed products at the discretion of WillowWood. If a replacement product is provided to a customer, the warranty period will be the number of days remaining on the original warranty. WillowWood's liability shall not exceed the purchase price of the product. WillowWood is not responsible for any indirect, incidental or consequential damage.

Retention of Rights

WillowWood retains all intellectual property rights reflected in its physical products, regardless of the transfer of the physical products to another party or parties.

Trial Period

The trial period for an Alpha Liner is 30 days from the date of the seller's invoice. The policy for Alpha Liners returned during the trial period is as follows:

- Credit will be offered for patient rejection, defined as "product not appropriate for patient"
- For all other circumstances, exchange for a different Alpha Liner will be offered

Returns

A customer's preferred distributor will handle all liner returns.

If you wish to return a product purchased directly from WillowWood, follow these simple steps:

A. Be prepared to provide the following information:

- Account # or company name
- Purchase order # and/or serial #
- Quantity to be returned
- Reason for return

B. Please call Customer Care at 1.800.848.4930 to get a Return Merchandise Authorization (RMA) number. This must be done before the product is returned.

C. When packing the product, place the RMA number prominently on the outside and inside of the package containing the product being returned.

By following these steps you will guarantee that your return will be processed expeditiously and credited accurately.

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WillowWood Global LLC
Ohio Willow Wood Company B.V.