



META® Foot Shell

Product Description

The META®Foot Shell is a cosmetic and protective component designed for use with compatible WillowWood META prosthetic feet that provide fit in a variety of footwear. The META Foot Shell is intended for use by individuals with lowerlimb amputations who have been prescribed a compatible WillowWood prosthetic foot.

Components Included

- META Foot Shell

Instructions

- The size and side of the Foot Shell and prosthetic foot must match
- Wear WillowWood prosthetic feet with a Foot Shell, Spectra Sock, and appropriate footwear. Please follow all instructions of the specific foot model.
- Only use a Foot Shell Tool for installation and removal of Foot Shell
- Secure prosthetic foot under heel tab of Foot Shell

Contraindications

- Use with noncompatible prosthetic feet.
- Use of a damaged foot shell.
- Any configuration not approved by the prescribing prosthetist.

Cleaning and Maintenance

- Inspect the foot shell regularly for cracks, excessive wear, or permanent deformation.
- Clean the foot shell using mild soap and water.
- Allow all components to air dry completely before reassembly.
- If exposed to chlorinated/salt water, solvents, or abrasives, immediately disassemble, rinse with water, and dry.
- Prolonged exposure to water or UV light may cause damage.
- Exposure to dyed clothing or footwear may cause discoloration

Important Safety Information

Do not continue use if the foot shell becomes damaged, loose, or excessively worn. Contact a certified prosthetist for evaluation and replacement if necessary.

WILLOWWOOD®

Warranty

The warranty for the Foot Shell is sixmonths from date of patient fitting. Use of the Foot Shell for amputees whose modified body weight is more than 166 kg (365 lbs) or who engage in extremely high and abusive activity is against WillowWood's recommendations and will void the sixmonth warranty. Modified body weight is defined as the weight of the amputee plus any loads carried by the amputee. "Extremely high and abusive activities" are defined as activities such as skydiving, karate, and judo; activities that could result in injury to an individual's natural feet; and activities that expose the prosthesis to corrosives such as salt water without adequate rinsing following exposure.

Warranty Disclaimer

WillowWood warrants that each product manufactured will, at the time of delivery, be of workmanlike quality and substantially free of defects. WILLOWWOOD MAKES NO OTHER WARRANTY, IMPLIED, OR EXPRESSED, AND MAKES NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. This warranty shall terminate immediately upon an action to combine our products with other materials or in any manner to change the nature of our products. The sole remedy is replacement of the products or credit for the products. WillowWood's liability shall not exceed the purchase price of the product. WillowWood shall not be liable for any indirect, incidental, or consequential damage.

WillowWood Retention of Rights

WillowWood retains all intellectual property rights reflected or incorporated in its physical products, regardless of the transfer of the physical products to another party or parties.

Reporting a Serious Incident

In the event of a serious incident occurring as a direct consequence of using this medical device, which either directly or indirectly caused (a) the death of the patient/user or another person, (b) the temporary or permanent serious deterioration of the patient's, user's, or other person's state of health, or (c) a serious public health threat, please immediately report the incident by contacting WillowWood Global or, if in the EU, the authorized European Union representative and the competent authority within the member state in which the patient/user is established.

WILLOWWOOD®

WillowWood Global LLC
15441 Scioto Darby Road
Mt. Sterling, OH 43143
téléphone 740.869.3377 / 800.848.4930
télécopieur 740.869.4374 www.willowwood.com
Brevet www.willowwood.com/patent



Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands



Electronic Instructions for Use are available at www.willowwood.com
A paper copy of the Instructions for Use is available upon request at
no additional cost.

