


Ordering Guidelines: Alpha[®] Silicone Liners



For liners with a progressive profile, the size of the accordion distal umbrella changes with the size of the liner.

SIZING GUIDE

Please refer to the illustrations for recommendations for fitting Alpha Silicone Liners.

TIP


If the limb circumferences fall into the range of both the Medium Plus 2 and the Large 1:

- choose **Medium Plus 2** if the distal end is bony and tapered
- choose **Large 1** if the distal end is broad and cylindrical

	Medium 1	Medium 2	Medium Plus 1	Medium Plus 2	Large 1	Large 2	Large Plus	Extra Large
Proximal Measurement	26 - 44 cm (10 - 17 in)	29 - 44 cm (11.5 - 17 in)	30 - 54 cm (12 - 21 in)	33 - 54 cm (13 - 21 in)	33 - 64 cm (13 - 25 in)	36 - 64 cm (14 - 25 in)	41 - 80 cm (16 - 31.5 in)	46 - 90 cm (18 - 35.5 in)
Distal Measurement	20 - 24 cm (8 - 9.5 in)	24 - 29 cm (9.5 - 11.5 in)	20 - 24 cm (8 - 9.5 in)	24 - 29 cm (9.5 - 11.5 in)	28 - 32 cm (11 - 12.5 in)	32 - 36 cm (12.5 - 14 in)	32 - 39 cm (12.5 - 15.5 in)	37 - 45 cm (14.5 - 18 in)

RETROFITTING

Patients wearing Alpha Hybrid Liners with a progressive profile may be retrofit with Alpha Silicone Liners with a progressive profile in the same base size (Medium Plus to Medium Plus, or Large to Large).



*A Medium Plus 3 size is available upon request to retrofit patients that are near the upper limit of the Alpha Hybrid Liner Medium Plus size range.

ITEM NUMBERS

Size	Locking Item #	Cushion Item #
Medium 1	S352-6590-1	S350-6590-1
Medium 2	S352-6590-2	S350-6590-2
Medium Plus 1	S352-6594-1	S350-6594-1
Medium Plus 2	S352-6594-2	S350-6594-2
Medium Plus 3	S352-6594-3	S350-6594-3
Large 1	S352-6596-1	S350-6596-1
Large 2	S352-6596-2	S350-6596-2
Large Plus	S352-6598-1	S350-6598-1
Extra Large	S352-6597-1	S350-6597-1

L-CODE SUGGESTIONS

L5679 Cushion Liner (bill once per liner)
L5673 Locking Liner (bill once per liner)

WARRANTY

12 months from date of invoice

TRIAL PERIOD

30 days from date of invoice