

*WillowWood®*

# ALPHA CONTROL™ LINER SYSTEM

*A partnership between WillowWood and Coapt*



**Get more out of your  
Pattern Recognition System.**

---

The next-generation prosthetic interface solution for upper extremity pattern recognition users.

# ALPHA CONTROL™ LINER SYSTEM

## The Next-Generation Prosthetic Interface Solution for Upper Extremity Pattern Recognition Users

The Alpha Control Liner System was developed in collaboration with Coapt and Shirley Ryan AbilityLab and is designed exclusively for use with the Coapt Gen2 Complete Control System. The Alpha Control Liner provides an improved solution over standard myoelectric interface methods that simplifies the fabrication process and delivers an improved patient experience.



### EASE OF SITE SELECTION AND SOCKET FITTING

Practitioners experience a simplified process for selecting myo-sites on the residual limb and modifying a positive model – reducing both the number of fittings and adjustments required.


### EXPANDED FUNCTIONAL ENVELOPE

By ensuring consistent electrode contact and enhancing interface efficiency, the Alpha Control Liner System allows amputees to achieve greater control through an expanded movement range.


### PATIENT COMFORT AND EASE OF USE

The Alpha Control Liner System provides distributed suspension via the liner and lock – leading to comfortable socket designs that are easier to don and doff than existing myoelectric interface methods.

The Alpha Control Liner is available in three sizes (Extra Small, Small and Medium), with nine potential lead lengths (ranging from 6 cm to 22 cm, in 2 cm increments) – enabling simplified myo-site selection and accommodating the vast majority of transradial and transhumeral limb sizes.



## WHAT'S INCLUDED



### ALPHA CONTROL LINER SYSTEM

- Alpha Control Liner
- Fiber Optic Cable and LED Cap
- Electrodes
- Module
- Lock
- Lock Dummy

## HOW TO ORDER

### Alpha Control Liner System

Ordering example:  
Medium Liner with 8 cm lead length

# ACLS-049-M08

### Replacement Alpha Control Liner Only

Ordering example:  
Medium Liner with 8 cm lead length

# ACL-049-M08

Model	Liner Size	Lead Length in cm		Model	Liner Size	Lead Length in cm	
ACLS-049 = Alpha Control Liner System	X = Extra Small S = Small M = Medium	06 cm 08 cm 10 cm 12 cm 14 cm	16 cm 18 cm 20 cm 22 cm	ACL-049 = Alpha Control Liner-Replace- ment Liner	X = Extra Small S = Small M = Medium	06 cm 08 cm 10 cm 12 cm 14 cm	16 cm 18 cm 20 cm 22 cm

### ALPHA CONTROL SYSTEM PART NUMBERS

Leads	Extra Small	Small	Medium
6 cm	ACLS-049-X06	ACLS-049-S06	ACLS-049-M06
8 cm	ACLS-049-X08	ACLS-049-S08	ACLS-049-M08
10 cm	ACLS-049-X10	ACLS-049-S10	ACLS-049-M10
12 cm	ACLS-049-X12	ACLS-049-S12	ACLS-049-M12
14 cm	ACLS-049-X14	ACLS-049-S14	ACLS-049-M14
16 cm	ACLS-049-X16	ACLS-049-S16	ACLS-049-M16
18 cm	ACLS-049-X18	ACLS-049-S18	ACLS-049-M18
20 cm	ACLS-049-X20	ACLS-049-S20	ACLS-049-M20
22 cm	ACLS-049-X22	ACLS-049-S22	ACLS-049-M22

### ALPHA CONTROL LINER PART NUMBERS

Leads	Extra Small	Small	Medium
6 cm	ACL-049-X06	ACL-049-S06	ACL-049-M06
8 cm	ACL-049-X08	ACL-049-S08	ACL-049-M08
10 cm	ACL-049-X10	ACL-049-S10	ACL-049-M10
12 cm	ACL-049-X12	ACL-049-S12	ACL-049-M12
14 cm	ACL-049-X14	ACL-049-S14	ACL-049-M14
16 cm	ACL-049-X16	ACL-049-S16	ACL-049-M16
18 cm	ACL-049-X18	ACL-049-S18	ACL-049-M18
20 cm	ACL-049-X20	ACL-049-S20	ACL-049-M20
22 cm	ACL-049-X22	ACL-049-S22	ACL-049-M22

*\*Liner is only available in a 3 mm uniform profile*

### ALPHA CONTROL LINER ACCESSORY PART NUMBERS

Part Number	Description
ACL-300	Alpha Control Liner Electrode Pack of 13
ACL-301	Alpha Control Liner Lock (with Button)
ACL-302	Alpha Control Module
ACL-303	Alpha Control Liner Cap
ACL-304	Alpha Control Liner Fabrication Dummy Kit
ACL-305	Alpha Control Lock Button (Replacement)

**FOR MORE INFORMATION  
CALL 800.848.4930 OR 740.869.3377  
OR VISIT WILLOWWOOD.COM**

