


Alpha® Silicone Liners Sizing Guide

Please refer to the illustrations for recommendations for fitting Alpha Silicone Liners. Shown below are limb circumferences that the various liners will fit, measured at 4 cm (1.5 in) and 30 cm (12 in) from the distal end.

If converting from Alpha Classic or Alpha Hybrid Liners with a progressive profile, use the two measurements to determine the appropriate size in Alpha Silicone Liners. Patients may be retrofit with Alpha Silicone Liners with a progressive profile in the same base size (Medium Plus to Medium Plus, or Large to Large).



NOTE: Accordion distal umbrella sizes change with each size liner.

	Medium 1	Medium 2	Medium Plus 1	Medium Plus 2	Large 1	Large 2	Large Plus	Extra Large
Proximal Measurement	26 - 44 cm (10 - 17 in)	29 - 44 cm (11.5 - 17 in)	30 - 54 cm (12 - 21 in)	33 - 54 cm (13 - 21 in)	33 - 64 cm (13 - 25 in)	36 - 64 cm (14 - 25 in)	41 - 80 cm (16 - 31.5 in)	46 - 90 cm (18 - 35.5 in)
Distal Measurement	20 - 24 cm (8 - 9.5 in)	24 - 29 cm (9.5 - 11.5 in)	20 - 24 cm (8 - 9.5 in)	24 - 29 cm (9.5 - 11.5 in)	28 - 32 cm (11 - 12.5 in)	32 - 36 cm (12.5 - 14 in)	32 - 39 cm (12.5 - 15.5 in)	37 - 45 cm (14.5 - 18 in)
Proximal Measurement 30 cm (12 in) from distal								
Distal Measurement 4 cm (1.5 in) from distal								
	Medium 1	Medium 2	Medium Plus 1	Medium Plus 2	Large 1	Large 2	Large Plus	Extra Large
Part Numbers								
Locking	S352-6590-1	S352-6590-2	S352-6594-1	S352-6594-2	S352-6596-1	S352-6596-2	S352-6598-1	S352-6597-1
Cushion	S350-6590-1	S350-6590-2	S350-6594-1	S350-6594-2	S350-6596-1	S350-6596-2	S350-6598-1	S350-6597-1