Alpha SmartTemp™ AK Liner Sizing Guide

Please refer to the illustrations for recommendations for fitting Alpha SmartTemp AK Liners. Shown below are limb circumferences that the various liners will fit, measured at 4 cm (1.5 in) and 30 cm (12 in) from the distal end.

NOTE: Accordion distal umbrella sizes change with each size liner.

Proximal Measurement	30 - 54 cm (12 - 21 in) 20 - 24 cm	33 - 54 cm (13 - 21 in)	33 - 64 cm (13 - 25 in)	36 - 64 cm (14 - 25 in)	41 - 80 cm (16 - 31.5 in)	46 - 90 cm (18 - 35.5 in)
			(13 - 25 in)	(14 - 25 in)	(16 - 31.5 in)	(18 - 35 5 in)
	20 - 24 cm					(10 33.3 111)
Distal Measurement	(8 - 9.5 in)	24 - 29 cm (9.5 - 11.5 in)	28 - 32 cm (11 - 12.5 in)	32 - 36 cm (12.5 - 14 in)	32 - 39 cm (12.5 - 15.5 in)	37 - 45 cm (14.5 - 18 in)
Proximal Measurement 30 cm (12 in) from distal	M ₊₁	M ₊₂		L ₂		
Distal Measurement 4 cm (1.5 in) from distal						

