

Date: MM.DD.YYYY

**CUSTOMER INFORMATION**

Clinician Name: \_\_\_\_\_

Facility/Clinic: \_\_\_\_\_

Preferred Contact Method:

Bill to Account #: \_\_\_\_\_ PO #: \_\_\_\_\_

Phone: \_\_\_\_\_

Shipping Address Line 1: \_\_\_\_\_

Email Address: \_\_\_\_\_

Shipping Address Line 2: \_\_\_\_\_

**SHIPPING INFORMATION**

Carrier: Method:

- UPS     Ground     3 - Day     2 - Day  
 FedEx     Next Day     Next Day Saver

Patient ID/Last Name: \_\_\_\_\_

Amputation Level: Amputation Side:

- TT     Left     Bilateral  
 Symes     Right

\* Requires separate order forms. Indicate side on each form.

**DUPLICATE LINER INFORMATION**


Duplicate Liner Serial Number: \_\_\_\_\_ Qty: \_\_\_\_\_

Cushion - ALC-DES-DP

Locking - ALL-DES-DP

DUO - ADC-DES-DP

**NEW LINER INFORMATION** (Fill in for new liners or alterations to existing liners)



Measurements:  
 cm     in

**LINER QUANTITY**

Qty: \_\_\_\_\_

Additional Liner codes:  
Cushion: ALC-DES-DP  
Locking: ALL-DES-DP  
DUO: ADC-DES-DP

**LINER TYPE**  
(select one)

- Cushion  
ALC-DES-EO  
 Classic  
 Hybrid

- Locking  
ALL-DES-EO  
 Classic  
 Hybrid

- Duo  
ADC-DES-EO

**FABRIC TYPE**  
(select one)

- Original  
 Buff  
 Gray/Green  
 Brown

- Spirit  
 Buff  
 Gray/Green  
 Brown

- Max  
Buff  
 Select  
 Gray  
 Taupe

**LINER DESIGN**  
(select one)

Would you like a Design Technician to contact you before starting this order?  
 Yes     No

Requested Starting Thickness

\_\_\_\_\_ mm Base Liner Thickness  
4mm to 12mm ( 6mm Standard )

\_\_\_\_\_ mm Distal Thickness  
4mm to 25mm ( 9mm Standard, 4 mm for Symes )

\_\_\_\_\_ mm Proximal Thickness  
3mm to 9mm ( 4mm Standard )

\_\_\_\_\_ mm Medial Thickness  
3mm to 9mm ( 4mm Standard )

Fill Invagination  
 Yes (recommended)     No

Stovepipe Liner  
 Yes (bulbous distal)     No

**SHAPE CREATION:**  
(select one)

- Cast  
 STL / AOP  
 Omega / Tracer  
 CAD File Name: \_\_\_\_\_  
 Cast Shipped Separately

Notes:

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_