

Date: MM.DD.YYYY

CUSTOMER INFORMATION

Clinician Name: _____
 Preferred Contact Method:
 Phone: _____
 Email Address: _____

Facility/Clinic: _____
 Bill to Account #: _____ PO #: _____
 Shipping Address Line 1: _____
 City/State/Zip Code _____

SHIPPING INFORMATION

Carrier: _____ Method: _____
 UPS Ground 3 - Day 2 - Day
 FedEx Next Day Next Day Saver

PATIENT INFORMATION


Patient ID/Last Name: _____
 Amputation Level: _____ Amputation Side:
 TF Left Bilateral
 KD Right
* Requires separate order forms. Indicate side on each form.

DUPLICATE LINER INFORMATION

Duplicate Liner Serial Number: _____ Qty: _____
Cushion - ALC-DES-DP Locking - ALL-DES-DP DUO - ADC-DES-DP

Starting July 1st 2024, WW will charge a \$100 service fee for Design Assessment and New Custom Model Creation - for any duplicate Alpha DESIGN liner orders placed separate from an initial (EO or OT DESIGN Liner Order). If a duplicate is ordered simultaneously with an initial Alpha DESIGN Liner order, the service fee will not be levied

NEW LINER INFORMATION (Fill in for new liners or alterations to existing liners)



Measurements: cm in

LINER TYPE:
(select one)

FABRIC TYPE
(select one)

LINER DESIGN
(select one)

- Cushion
ALC-DES-EO
- Classic
- Hybrid

- Original
Buff
- Spirit
Buff
- Max
Buff
- Select
 - Gray
 - Taupe

Would you like a Design Technician to contact you before starting this order?
 Yes No

LINER QUANTITY

Qty: _____

Additional Liner codes:
Cushion: ALC-DES-DP
Locking: ALL-DES-DP
DUO: ADC-DES-DP

- Locking
ALL-DES-EO
- Classic
- Hybrid

- Duo
ADC-DES-EO

SHAPE CREATION:
(select one)

- Cast
- STL / AOP
- Omega / Tracer
- CAD File Name: _____
- Cast Shipped Separately

Requested Starting Thickness

_____ mm Base Liner Thickness
4mm to 12mm (6mm Standard)

_____ mm Distal Thickness
4mm to 25mm (9mm Standard, 4 mm for Symes)

_____ mm Proximal Thickness
3mm to 9mm (4mm Standard)

_____ mm Medial Thickness
3mm to 9mm (4mm Standard)

Fill Invagination
 Yes (recommended) No

Stovepipe Liner
 Yes (bulbous distal) No

Notes:

