



THE D.I.Y. Adjustable, Thermoformable, Definitive Socket

Why XtremityTT?

HEAT, SHAPE, ASSEMBLE

Fabricate a definitive socket in 1 hour.

QUICK & EFFECTIVE ADJUSTMENTS

Make quick and effective post-delivery adjustments in minutes.

Simple Components



Socket Preform & Base Plate



Pin-Lock



Suction



Vacuum



Socket Advantages

- Adjustable Base Plate: 16mm slide, 30° rotation.
- Modular, Interchangeable Suspension: Pin-lock, Suction, and Vacuum.
- Improved Patient Comfort: Enhanced material properties for flexibility.*
- Tested: Socket and componentry tested to ISO 10328 Standards.*
- Seven sizes available
- 275lb weight limit



Health & Safety

Socket fabrication prioritizes health and safety, eliminating the carbon itch when trimming and grinding.

* References available upon request.

Contact your WW Account Manager

The D.I.Y.,
Flexible Inner Socket
with a consistent wall
thickness throughout



Eliminates grinding
and remakes due to
thickness errors.



Fabricates in less
than one hour.



Minimal lab experience
required: Single person
fabrication with
consistent results.



Heating in the same
Benchtop Heating Unit
as XtremityTT socket.

Scan here for
fabrication videos.

